

ZNF131 Rabbit pAb

货号: B18502

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	67kDa/71kDa
实测分子量	80-85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse kidney,Mouse lung,Rat testis,Rat kidney
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-240 of human ZNF131 (NP_003423.1).
序列	MAFRHLIEFTYTAKLMIQGEEEEANDVWKAAEFLQMLEAIKALEVRNKEN SAPLEENTTGKNEAKKRKIAETSNVITESLP SAESEPVEIEVEIAEGTIEVEDEGIETLEEVASAKQSVKYIQSTGSSDD SALLADITSKYRQGD RKGQIKEDGCPSDPT

靶点信息

研究背景	Plays a role during development and organogenesis as well as in the function of the adult central nervous system (By similarity. May be involved in transcriptional regulation as a repressor of ESR1/ER-alpha signaling.
基因ID	7690

基因名	ZNF131
Swiss	P52739
别名	ZNF131;ZBTB35;pHZ-10

产品验证



Western blot analysis of ZNF131 expressed in Mouse kidney,Rat testis,HeLa using ZNF131 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn