

ZNF19 Rabbit pAb

货号: B18506

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	51kDa/52kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ZNF19 (NP_008892.2).
序列	MAAMPLKAQYQEMVTFEDVAVHFTKTEWTGLSPAQRALYRSVMLENFGNLTALGYPVPKPALISLLERGDMAWGLEAQD DPPAERTKNVCKDVETNIDSESTLIQGISEERDGMMSHGQLKSVQRTDFPETRNVEKHQDIPTVKNIQGKVPRIPCARKP FICEECGKSFSYFSYYARHQ

靶点信息

研究背景	The protein encoded by this gene contains a zinc finger, a nucleic acid-binding domain present in many tr anscription factors. This gene is located in a region next to ZNF23, a gene also encoding a zinc finger prot ein, on chromosome 16.
基因ID	7567

基因名	ZNF19
Swiss	P17023
别名	ZNF19;KOX12

产品验证



Western blot analysis of ZNF19 expressed in U-251MG using ZNF19 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn