

# KCNJ15 Rabbit pAb

货号: B18539

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IHC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:200 - 1:2000<br><b>IHC:</b> 1:50 - 1:200                                                |
| 理论分子量 | 42kDa                                                                                                |
| 实测分子量 | 35kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | 293T,Mouse liver,Mouse kidney,Rat liver                                                              |
| 细胞定位  | Membrane,Multi-pass membrane protein                                                                 |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                              |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 221-375 of human KCNJ 15 (NP_002234.2).                                        |
| 序列   | VTKEGERILLNQATVKFHVDDSSSESPFLILPMTFYHVLDETSPLRDLTPQNLKEKEFELVLLNATVESTSAVCQSRTSYIPEEI YWGFEFVPVVSLSKNGKYVADFSQFEQIRKSPDCTFYCADSEKQQLLEKYRQEDQRERELRTLLLQGSNV |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein has a greater tendency to allow potassium to flow into a cell rather than out of a cell. Eight transcript variants encoding the same protein have been found for this gene. |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                           |
|-------|---------------------------|
| 基因ID  | 3772                      |
| 基因名   | KCNJ15                    |
| Swiss | Q99712                    |
| 别名    | KCNJ15;IRKK;KIR1.3;KIR4.2 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)