

# SIK1 Rabbit pAb

货号: B18563

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                            |
| 理论分子量 | 84kDa                                                                                                |
| 实测分子量 | 80kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Jurkat,SGC-7901,Mouse stomach                                                                        |
| 细胞定位  | Cytoplasm,Nucleus                                                                                    |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                         |
|------|---------------------------------------------------------------------------------------------------------|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SIK1 (NP_775490.2).   |
| 序列   | MVIMSEFSADPAGQGQGGQKPLRVGFYDIERTLGKGNFAVVKLARHRVTKTQVAIKIIDKTRLDSSNLEKIYREVQLMKLLN<br>HPHIKLYQVMETKDMLY |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a serine/threonine protein kinase that contains a ubiquitin-associated (UBA) domain. The encoded protein is a member of the adenosine monophosphate-activated kinase (AMPK) subfamily of kinases that play a role in conserved signal transduction pathways. A mutation in this gene is associated with early infantile epileptic encephalopathy 30. |
| 基因ID | 150094                                                                                                                                                                                                                                                                                                                                                                 |

|       |                                  |
|-------|----------------------------------|
| 基因名   |                                  |
| Swiss | P57059                           |
| 别名    | SIK1;MSK;SIK;SIK-1;SNF1LK;EIEE30 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)