

FARS2 Rabbit pAb

货号: B18566

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	52kDa
实测分子量	52kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A549,rat kidney
细胞定位	Mitochondrion matrix
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human FARS2 (NP_006558.1).
序列	MVGSALRRGAHAYVYLVS KASHISRGHQHQAWGSRPPAAECATQ RAPGSVVELLGKSYPQDDHSNLTRKVLTRVGRNL HNQQHHPLWLIKERVKEHFYKQYVGRFGTPLFSVYDNLSPVVTTWQNFDSLIPADHPSRKKGDNYLNRTHMLRAHTSA HQWDL LHAGLDAFLVVG DVYRRDQIDSQHYP IFHQLEAVRLFSKHEL FAGIKDGE SLQLFEQSSRSAHKQETH TMEAVKL VEFDLKQTLTRL

靶点信息

研究背景	This gene encodes a protein that transfers phenylalanine to its cognate tRNA. This protein localizes to the mitochondrion and plays a role in mitochondrial protein translation. Mutations in this gene can cause combined oxidative phosphorylation deficiency 14 (Alpers encephalopathy). Alternative splicing results in multiple transcript variants.
------	---

基因ID	10667
基因名	FARS2
Swiss	O95363
别名	FARS2;COXPD14;FARS1;HSPC320;PheRS;SPG77

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn