

TEKT1 Rabbit pAb

货号: B18569

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	Cytoplasm,cilium axoneme,cytoskeleton,flagellum axoneme
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TEKT1 (NP_444515.1).
序列	MAKLLQPPPKFLPSEWHIANKNQYHRADAQRSRSERLVAESQRLVDEIEKTTRKSQSDVNKKLEQRLEEVQFWKKELDDKLEQLVNVTDLLIYKIRLEKALETLKEPLHITETCLAYREKRIGIDLVDHDTVEHELIEAEIIQGIMALLTRTLEEASEQIRMNRSAYNLEKDLKDKFVALTIDDICFSLNNNSPNIRYSENAVRIEPNSVSLDWLDFSSSTNVEKADKQRNNSMLKALVDRILSQTANDLRKQCDVVDTAfKNGlKDKDARDKLADHLAKVMEEIASQEKNTALE

靶点信息

研究背景	This gene product belongs to the tektin family of proteins. Tektins comprise a family of filament-forming proteins that are coassembled with tubulins to form ciliary and flagellar microtubules. This gene is predominantly expressed in the testis and in mouse, tektin 1 mRNA was localized to the spermatocytes and round spermatids in the seminiferous tubules, indicating that it may play a role in spermatogenesis.
------	--

基因ID	83659
基因名	TEKT1
Swiss	Q969V4
别名	TEKT1;tektin-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn