

VEPH1 Rabbit pAb

货号: B18572

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	19kDa/23kDa/89kDa/94kDa
实测分子量	95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,Mouse kidney
细胞定位	Cell membrane,Cytoplasmic side,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 448-547 of human VEPH 1 (NP_078897.2).
序列	KSLAFHTMLTKGVGSDDGEDENRGDIPASISLSEIDPLGQGNDKLPFKTDTERSQLGESSVSYPNIIHIDSENLSETVKENS QEETPETTASPIEYQDKL

靶点信息

研究背景	Interacts with TGF-beta receptor type-1 (TGFB R1, impeding its nuclear accumulation and resulting in impaired TGF-beta signaling. May also affect FOX O, Hippo and Wnt signaling.
基因ID	79674

基因名	VEPH1
Swiss	Q14D04
别名	VEPH1;MELT;VEPH

产品验证



Western blot analysis of VEPH1 expressed in A-431,Mouse kidney using VEPH1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn