

HSPA12B Rabbit pAb

货号: B18583

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	75kDa
实测分子量	72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 552-686 of human HSPA12B (NP_443202.3).
序列	VPGRHPPEKLLVRDGRRWCTDVFERFVAAEQSVALGEEVRRSYCPARPGQRRVLINLYCCAEDARFITDPGVRKCGALSL ELEPADCGQDTAGAPPGRREIRAAMQFGDTEIKVTAVDVSTNRSVRASIDFLSN

靶点信息

研究背景	The protein encoded by this gene contains an atypical heat shock protein 70 (Hsp70) ATPase domain and is therefore a distant member of the mammalian Hsp70 family. This gene may be involved in susceptibility to atherosclerosis. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	116835

基因名	HSPA12B
Swiss	Q96MM6
别名	HSPA12B;C20orf60

产品验证



Western blot analysis of HSPA12B expressed in LO2 using HSPA12B Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn