

EIF4ENIF1 Rabbit pAb

货号: B18597

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	88kDa/108kDa
实测分子量	178kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,HT-29,LO2
细胞定位	Cytoplasm,Nucleus,Nucleus speckle,PML body
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human EIF4ENIF1 (NP_001157973.1).
序列	LGRRRIGSGRIISARTFEKDHRLSDKDLRDLRDRERDFKDKRFRREFGDSKRVFGERRRRNDSYTEEEPEWFSAGPTSQSETIELTGFDKILEEDHKGR

靶点信息

研究背景	The protein encoded by this gene is a nucleocytoplasmic shuttle protein for the translation initiation factor eIF4E. This shuttle protein interacts with the importin alpha-beta complex to mediate nuclear import of eIF4E. It is predominantly cytoplasmic; its own nuclear import is regulated by a nuclear localization signal and nuclear export signals. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	56478
基因名	EIF4ENIF1
Swiss	Q9NRA8
别名	EIF4ENIF1;4E-T;Clast4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn