

STRBP Rabbit pAb

货号: B18599

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	71kDa/73kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 40-190 of human STRBP (NP_060857.2).
序列	CALKHVSDWLDETNGTKTEGETEVKKDEAGENYSKDQGGRTLCGVMRIGLVAKGLLIKDDMDLELVLMMCKDKPTETLLNTVKDNLPIQIQKLTEEKYQVEQCVNEASIIIRNTKEPTLTLKVILTSPLIRDELEKKDGENVSMKDPPDLL

靶点信息

研究背景	Involved in spermatogenesis and sperm function. Plays a role in regulation of cell growth. Binds to double-stranded DNA and RNA. Binds most efficiently to poly(I:C RNA than to poly(dI:dC DNA. Binds also to single-stranded poly(G RNA. Binds non-specifically to the mRNA PRM1 3'-UTR and adenovirus VA RNA (By similarity).
基因ID	55342

基因名	STRBP
Swiss	Q96SI9
别名	STRBP;HEL162;ILF3L;SPNR;p74

产品验证



Western blot analysis of STRBP expressed in HeLa using STRBP Rabbit pAb at 1:1000. Secondary antibody : HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn