

FOXB1 Rabbit pAb

货号: B18608

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	34kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human FOX B1 (NP_036314.2).
序列	LPNQLTTMGSSSLGTGWPHVYGSAGMIDSATPISMASGDYSAYGVPLKPLCHAAGQTLPAIPVPIKPTPAAVPALPALPAPIPTLLSNSPPSLSPTSSQTATSQSSPATPSETLTSPASALHSVAVH

靶点信息

研究背景	Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in several processes, including mammary gland development; nervous system development; and visual learning. Predicted to act upstream of or within spinal cord development; thalamus development; and urogenital system development. Predicted to be located in nucleus. Predicted to be part of chromatin.
基因ID	27023

基因名	FOXB1
Swiss	Q99853
别名	FOXB1;FKH5;HFKH-5

产品验证



Western blot analysis of FOXB1 expressed in HepG2 using FOXB1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn