

TRPC6 Rabbit pAb

货号: B18628

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	93kDa/100kDa/106kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SGC996,Mouse lung,Mouse brain,Rat brain
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 732-931 of human TRPC 6 (NP_004612.2).
序列	QEIEDDADVEWKFAKALWFSYFEEGRTPVPFNLVPSPKSLFYLLLKLKKWISELFQGHKKGFQEDAEMNKINEEKKLGIL GSHEDL SKLSLDKKQVGHNKQPSIRSSDFHLNSFNPPRQYQKIMKRLIKRYVLQAQIDKESDEVNEGELKEIKQDISSLR YELLEESQNTEDLAELIRELGEKLSMEPNQEETNR

靶点信息

研究背景	The protein encoded by this gene forms a receptor-activated calcium channel in the cell membrane. The channel is activated by diacylglycerol and is thought to be under the control of a phosphatidylinositol second messenger system. Activation of this channel occurs independently of protein kinase C and is not triggered by low levels of intracellular calcium. Defects in this gene are a cause of focal segmental glomerulosclerosis 2 (FSGS2).
------	---

基因ID	7225
基因名	TRPC6
Swiss	Q9Y210
别名	TRPC6;FSGS2;TRP6

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn