

SIM2 Rabbit pAb

货号: B18645

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	63kDa/73kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HepG2
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 342-525 of human SIM2 (NP_033664.2).
序列	LSLEQVSTAKSQDSWRTALSTSQETRKLVKPKNTKMKTCLRTPYPPQQYSSFQMDKLECGQLGNWRASPPASAAAPPE LQPHSESSDLLYTPSYSLPFSYHYGHFPLDSHVFSKKPMPLPAKFGQPQGSPCEVARFFLSTLPASGECQWHYANPLVPSS SSPAKNPPEPPANTARHSLVPSYE

靶点信息

研究背景	This gene represents a homolog of the Drosophila single-minded (sim) gene, which encodes a transcription factor that is a master regulator of neurogenesis. The encoded protein is ubiquitinated by RING-IBR-RING-type E3 ubiquitin ligases, including the parkin RBR E3 ubiquitin protein ligase. This gene maps within the so-called Down syndrome chromosomal region, and is thus thought to contribute to some specific Down syndrome phenotypes. Alternative splicing of this gene results in multiple transcript variants.
------	--

基因ID	6493
基因名	SIM2
Swiss	Q14190
别名	SIM2;HMC13F06;HMC29C01;SIM;bHLHe15

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn