

# DAPK3 Rabbit pAb

货号: B18664

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IHC IF/ICC                                                                                        |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:100<br><b>IF/ICC:</b> 1:50 - 1:100                 |
| 理论分子量 | 37kDa/52kDa                                                                                          |
| 实测分子量 | 55kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | mouse brain,mouse lung,rat brain                                                                     |
| 细胞定位  | Cytoplasm,Cytoplasm,Nucleus                                                                          |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                          |
|------|----------------------------------------------------------------------------------------------------------|
| 抗原信息 | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DAPK3 (NP_001339.1).   |
| 序列   | MSTFRQEDVEDHYEMGEELGSGQFAIVRKCRQKGTGKEYAAKFIKKRRLSSSRRGVSREEIEREVNILREIRHPNIITLHDIFE<br>NKTDVVLILELVSGGE |

靶点信息

|      |                                                                                                                                                                                                             |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Death-associated protein kinase 3 (DAPK3) induces morphological changes in apoptosis when overexpressed in mammalian cells. These results suggest that DAPK3 may play a role in the induction of apoptosis. |
| 基因ID | 1613                                                                                                                                                                                                        |

|       |                    |
|-------|--------------------|
| 基因名   | DAPK3              |
| Swiss | O43293             |
| 别名    | DAPK3;DLK;ZIP;ZIPK |

产品验证



Western blot analysis of DAPK3 expressed in mouse brain,mouse lung,rat brain using DAPK3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)