

ARCN1 Rabbit pAb

货号: B18684

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	47kDa/57kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,HeLa,NIH/3T3,Mouse brain
细胞定位	COPI-coated vesicle membrane,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Golgi apparatus membra ne,Peripheral membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 212-511 of human ARC N1 (NP_001646.2).
序列	VAPAPARPSGSPKALKLGAKGKEVDNFVDKLGSEGETIMSSSMGKRTSEATKMHAPPINMESVHMKIEEKITLTCGRDGGGL QNMELHGMIMLRISDDKYGRIRLHVENEDKKGVQLQTHPNVDKKLFTAESLIGLKNPEKSFPVNSDVGVLKWRLQTTEES FIPLTINCWPSESGNGCDVNIEYELQEDNLELNDVVITIPLPSGVGAPVIGEIDGEYRHDSRRNTLEWCPLVIDAKNKSGSLE FSIAGQPNDFFPVQVSFVSKKNYCNIQVTKVTQVDGNSPVRFSTETTFLVDKYEIL

靶点信息

研究背景	This gene maps in a region, which include the mixed lineage leukemia and Friend leukemia virus integrati on 1 genes, where multiple disease-associated chromosome translocations occur. It is an intracellular pro tein. Archain sequences are well conserved among eukaryotes and this protein may play a fundamental r ole in eukaryotic cell biology. It has similarities to heat shock proteins and clathrin-associated proteins, a nd may be involved in vesicle structure or trafficking.
基因ID	372
基因名	ARCN1
Swiss	P48444
别名	ARCN1;COPD;SRMMD;archain 1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn