

B3GNT6 Rabbit pAb

货号: B18690

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa/42kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Mouse stomach,Mouse pancreas
细胞定位	Golgi apparatus membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human B3GN T6 (NP_619651.3).
序列	FLQAPRSPREERSPQEETPEGPTDAPAADEPPSELVPGPPCVANASANATADFEQLPARIQDFLRYRHCRHFLLWDAPAK CAGGRGVFLLAVKSAPEHYERRELIRRTWQGERSYGGRPVRRLFLLGTPGPEDEARAERLAELVALEAREHGDVLQWAF ADTFLNLTLK

靶点信息

研究背景	The protein encoded by this gene is a beta-1,3-N-acetylglucosaminyltransferase that adds an N-acetylglucosamine moiety to N-acetylgalactosamine-modified serine or threonine. The encoded enzyme is responsible for creating the core 3 structure of O-glycans, which are important components of mucin-type glycoproteins.
基因ID	192134

基因名	B3GNT6
Swiss	Q6ZMB0
别名	B3GNT6;B3Gn-T6;BGnT-6;beta-1;3-Gn-T6

产品验证



Western blot analysis of B3GNT6 expressed in HT-29,Mouse stomach,Mouse pancreas using B3GNT6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time : 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn