

DPPA2 Rabbit pAb

货号: B18699

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Jurkat,Mouse ovary,Mouse uterus,Rat ovary,Rat uterus
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human DPPA2 (NP_620170.3).
序列	MSDANLDSSKKNFLEGEVDDEESVILTLVPVKDDANMEQMEPSVSSTSDVKLEKPKKYNPGHLLQTNEQFTAPQKARCKI PALPLPTILPPINKVCRDTRLRDWCQQLGLSTNGKKIEVYLRRLHRHAYPEQRQDMPEMSQETRLQRCRKRKAVTKRARLQR SYEMNERAEETNTVEVITSAPGAMLASWARIAARAVQPK

靶点信息

研究背景	Binds to target gene promoters, including NKX2-5 and SYCE1, but not GATA4, and may be involved in the maintenance of the active epigenetic status of these genes.
基因ID	151871

基因名	DPPA2
Swiss	Q7Z7J5
别名	DPPA2;CT100;ECAT15-2;PESCRG1

产品验证



Western blot analysis of DPPA2 expressed in Mouse ovary,Rat ovary,SH-SY5Y using DPPA2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn