

MYNN Rabbit pAb

货号: B18726

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	10kDa/56kDa/65kDa/68kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3,A-549,HeLa,Mouse brain,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human MYNN (NP_061127.1).
序列	MQYSHHCEHLLERLNKQREAGFLCDCTIVIGEFQFKAHRNVLASFSEYFGAIYRSTSENNVFLDQSQVKADGFQKLEFIYT GTLNLDSWNVKEIHQAADYLKVEEVVTCKIKMEDFAFIANPSSTEISSITGNIELNQQTCLLTRDYNNREKSEVSTDLIQA NPKQGALAKKSSQTKKKKKAFNSPKTGQNKTQYPSDILENASVELFLDANKLPTPVVEQVAQINDNSEELTSVVENTFP AQDIVHTVTVKRKRKGSQPNCALKEHSMNSIASVKSPYEAENSGEELDQRYSKA

靶点信息

研究背景	This gene encodes a member of the BTB/POZ and zinc finger domain-containing protein family that are in volved in the control of gene expression. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 14.
基因ID	55892

基因名	MYNN
Swiss	Q9NPC7
别名	MYNN;OSZF;SBBIZ1;ZBTB31;ZNF902;myoneurin

产品验证



Western blot analysis of MYNN expressed in Mouse brain,OVCAR3,A-549 using MYNN Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn