

FGGY Rabbit pAb

货号: B18729

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	26kDa/27kDa/47kDa/50kDa/59kDa/62kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,A-431,Mouse kidney,Rat liver
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human FGGY (NP_060761.3).
序列	MSGGEQKPERYYVGVDVGTGSVRAALVDQSGVLLAFADQPIKNWEPQFNHHEQSSEDIWAACCVVTKKVQGIDLNQIRGLGFDATCSLVLDKQFHPLPVNQEGDSHRNVIMWLDHRAVSQVNRINETKHSVLQYVGGVMSVEMQAPKLLWLKENLREICWDKAGHFFDLPDFLSWKATGVTARSLCSLVCKWTYSAEKGWDDSFWKMIGLEDFVADNYSKIGNQVLP PGASLGNGLTPEAARDLGLLP

靶点信息

研究背景	This gene encodes a protein that phosphorylates carbohydrates such as ribulose, ribitol, and L-arabinitol. Genome-wide association studies in some populations have found an association between polymorphism s in this gene and sporadic amyotrophic lateral sclerosis, but studies of other populations have not been able to replicate this association. Alternative splicing results in multiple transcript variants.
------	---

基因ID	55277
基因名	FGGY
Swiss	Q96C11
别名	FGGY

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn