

GALNT10 Rabbit pAb

货号: B18730

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	31kDa/41kDa/61kDa/68kDa
实测分子量	69kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SGC-7901,OVCAR3,Mouse brain,Mouse testis,Rat heart,Rat brain
细胞定位	Golgi apparatus membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 444-603 of human GALNT10 (NP_938080.1).
序列	WDLPKFYPPVEPPAAAWGEIRNVGTGLCADTKHGALGSPLRLEGCVRGRGEAAWNNMQVFTFTWREDIRPGDPQHKKFCFDAISHTSPVTLYDCHSMKGNQLWKYRKDKTLYHPVSGSCMDCSESDHRIFMNTCNPSSLTQQWLFHTNSTVLEKFNRRN

靶点信息

研究背景	This gene encodes a member of the GalNAc polypeptide N-acetylgalactosaminyltransferases. These enzymes catalyze the first step in the synthesis of mucin-type oligosaccharides. These proteins transfer GalNAc from UDP-GalNAc to either serine or threonine residues of polypeptide acceptors. The protein encoded by this locus may have increased catalytic activity toward glycosylated peptides compared to activity toward non-glycosylated peptides.
------	---

基因ID	55568
基因名	GALNT10
Swiss	Q86SR1
别名	GALNT10;GALNACT10;PPGALNACT10;PPGANTASE10

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn