

Bif-1 Rabbit pAb

货号: B18741

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	29kDa/40kDa/43kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	Cytoplasm,Cytoplasmic vesicle,Golgi apparatus membrane,Midbody,Mitochondrion outer membrane,Peri pheral membrane protein,autophagosome membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 120-365 of human Bif-1 (NP_057093.1).
序列	KCGETQKRIGTADRELIQTSALNFLTPLRNFIEGDYKTIAKERKLLQNKRLDLDAAKTRLKKAKAAETRNSSEQELRITQSEF DRQAEITRLLLEGISSTHAHHLRCLNDFVEAQMTYYAQCYQYMLDLQKQLGSFPSNYLSNNNQTSVTPVPSVLPNAIGSSA MASTSGLVITSPSNLSDLKECSGSRKARVLYDYDAANSTELSLLADEVITVFSVVGMDSDWLMGERGNQKGKVPITYLELL N

靶点信息

研究背景	This gene encodes a SRC homology 3 domain-containing protein. The encoded protein interacts with the proapoptotic member of the Bcl-2 family, Bcl-2-associated X protein (Bax) and may be involved in regulat ing apoptotic signaling pathways. This protein may also be involved in maintaining mitochondrial morphol ogy. Alternate splicing results in multiple transcript variants.
------	---

基因ID	51100
基因名	SH3GLB1
Swiss	Q9Y371
别名	SH3GLB1;Bif-1;CGI-61;PPP1R70;dj612B15.2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn