

ERLIN1 Rabbit pAb

货号: B18769

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	38kDa
实测分子量	39kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,BxPC-3,DU145,OVCAR3,Mouse kidney,Rat testis
细胞定位	Endoplasmic reticulum membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 259-348 of human ERLIN1 (NP_006450.2).
序列	KADAEYYAAHKYATSNKHKLTPLEYLELKKYQAIASNSKIYFGSNIPNMFVDSSCALKYSDIRTGRESSLPKEALEPSGENVIQNKESTG

靶点信息

研究背景	The protein encoded by this gene is part of a protein complex that mediates degradation of inositol 1,4,5-trisphosphate receptors in the endoplasmic reticulum. The encoded protein also binds cholesterol and regulates the SREBP signaling pathway, which promotes cellular cholesterol homeostasis. Defects in this gene have been associated with spastic paraplegia 62.
------	--

基因ID	10613
基因名	ERLIN1
Swiss	O75477
别名	ERLIN1;C10orf69;Erlin-1;KE04;KEO4;SPFH1;SPG62;erlin-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn