

RBBP6 Rabbit pAb

货号: B18802

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	13kDa/106kDa/197kDa/201kDa
实测分子量	210KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562
细胞定位	Chromosome,Cytoplasm,Nucleus,centrosome,cytoskeleton,microtubule organizing center,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1240-1355 of human RBBP6 (NP_008841.2).
序列	SKVKQEKVKGKVRKVTGTGEGSSSTLVDTSTSSSTGGSPVRKSEKTDTKRTVIKTMEEYNNDNTAPAEDVIIMIQVPQSKWDKDDFESEEDVKSTQPISSVGKPAVKNVSTK

靶点信息

研究背景	The retinoblastoma tumor suppressor (pRB) protein binds with many other proteins. In various human cancers, pRB suppresses cellular proliferation and is inactivated. Cell cycle-dependent phosphorylation regulates the activity of pRB. This gene encodes a protein which binds to underphosphorylated but not phosphorylated pRB. Multiple alternatively spliced transcript variants that encode different isoforms have been found for this gene.
基因ID	5930

基因名	RBBP6
Swiss	Q7Z6E9
别名	RBBP6;MY038;P2P-R;PACT;RBQ-1;SNAMA

产品验证



Western blot analysis of RBBP6 expressed in K-562 using RBBP6 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn