

Smad6 Rabbit pAb

货号: B18815

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	26kDa/35kDa/53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,22Rv1,HeLa
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 257-496 of human Sma d6 (NP_005576.3).
序列	PTVCCNPYHFSRLCGPESPPPPYSRLSPRDEYKPLDLSSTLSYTETETATNSLITAPGEFSDASMSPDATKPSHWCSVAYW EHRTRVGRLYAVYDQAVSIFYDLPGSGFCLGQLNLEQRSESVRRTRSKIGFGILLSKEPDGVWAYNRGEHPFVNSPTLD APGGRALVVRKVPPGYSIKVFDFERSGLQHAPEPDAADGPYDPNSVRISFAKGWGPCYSRQFITSCPCWLEILLNNPR

靶点信息

研究背景	The protein encoded by this gene belongs to the SMAD family of proteins, which are related to Drosophila 'mothers against decapentaplegic' (Mad) and C. elegans Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein functions in the negative regulation of BMP and TGF-beta/activin-signalling. Multiple transcript variants have been found for this gene.
------	--

基因ID	4091
基因名	SMAD6
Swiss	O43541
别名	SMAD6;AOVD2;HsT17432;MADH6;MADH7

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn