

# CBFA2T3 Rabbit pAb

货号: B18835

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:100 - 1:500                                                                             |
| 理论分子量 | 7kDa/62kDa/63kDa/71kDa                                                                               |
| 实测分子量 | 75KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | K-562,Jurkat,Mouse spleen,Mouse thymus,Rat thymus                                                    |
| 细胞定位  | Golgi apparatus membrane,Nucleus,nucleolus,nucleoplasm                                               |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 305-390 of human CBFA 2T3 (NP_005178.4). |
| 序列   | DRDPLHPEHLSKR PCTLNPAQRYSPSNGPPQTPPPHYRLEDIAMAHHFRDAYRHPDPRELRRERHRPLVVPGSRQEEVIDH KLTER                               |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the myeloid translocation gene family which interact with DNA-bound tra nscription factors and recruit a range of corepressors to facilitate transcriptional repression. The t(16;21)( q24;q22) translocation is one of the less common karyotypic abnormalities in acute myeloid leukemia. Th e translocation produces a chimeric gene made up of the 5'-region of the runt-related transcription factor 1 gene fused to the 3'-region of this gene. This gene is also a putative breast tumor suppressor. Alternati ve splicing results in transcript variants. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                         |
|-------|-----------------------------------------|
| 基因ID  | 863                                     |
| 基因名   | CBFA2T3                                 |
| Swiss | O75081                                  |
| 别名    | CBFA2T3;ETO2;MTG16;MTGR2;RUNX1T3;ZMYND4 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)