

CPNE3 Rabbit pAb

货号: B18847

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	60kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	OVCAR3,Mouse brain,Mouse heart,Mouse kidney
细胞定位	Cell junction,Cell membrane,Cytoplasm,Nucleus,focal adhesion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 238-537 of human CPNE 3 (NP_003900.1).
序列	EASRSSPVEFECINEKKRQKKKSYKNSGVISVKQCEITVECTFLDYIMGGCQLNFTVGVDFTGSNGDPRSPDSLHYISPNG VNEYLTALWSVGLVIQDYDADKMFPAGFGAQIPQWQVSHEFPMNFNPSNPYCNGIQGIVEAYRSCLPQIKLYGPTNFSPIINHVARFAAAATQQQTASQYFVLLIITDGVITDLDETRQAIVNASRLPMSIIVGVGGADFSAMEFLDGDGGSLRSPLGEVAIRDIVQFVPPFRQFQNPAPKEALACVLAIEPQQVVGYPNTYKLLPPKNPATKQKQKQ

靶点信息

研究背景	Calcium-dependent membrane-binding proteins may regulate molecular events at the interface of the cell membrane and cytoplasm. This gene encodes a protein which contains two type II C2 domains in the amino-terminus and an A domain-like sequence in the carboxy-terminus. The A domain mediates interaction between integrins and extracellular ligands.
------	--

基因ID	8895
基因名	CPNE3
Swiss	O75131
别名	CPNE3;CPN3;PRO1071;copine-3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn