

CNKSR2 Rabbit pAb

货号: B18891

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	96kDa/101kDa/113kDa/117kDa
实测分子量	117kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Jurkat,HepG2
细胞定位	Cytoplasm,Membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 400-750 of human CNK SR2 (NP_055742.2).
序列	RKRFNIVEEDTVLYCYEYEKGRSSSQGRRESTPTYGKLRPISMPVEYNWVG DYEDPNKMKRDSRRENSLLRYMSNEKIAQE EYMFQRNSKKDTGKKSKKKGDKSNSPHYSLLPSLQMDALRQDIMGTPVPETTLYHTFQQSSLQHKS KKKKNGPIAGKSK RRISCKDLGRGDCEGWLWKKKDAKSYFSQKWKKYWFVLKDASLYWYINEEDEKAEGFISLPEFKIDRASECRKKYAFKAC HPKIKSFYFAAEHLDDMNRWLNRLNMLTAGYAERERIKQE QDYWSESDKEEADTPSTPKQDSPPPPYDTYPRPPSMSCAS PYVEAKHSRLSSTETSQSQSSHEEFRQEV T

靶点信息

研究背景	This gene encodes a multidomain protein that functions as a scaffold protein to mediate the mitogen-activated protein kinase pathways downstream from Ras. This gene product is induced by vitamin D and inhibits apoptosis in certain cancer cells. It may also play a role in ternary complex assembly of synaptic proteins at the postsynaptic membrane and coupling of signal transduction to membrane/cytoskeletal remodeling. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	22866
基因名	CNKSR2
Swiss	Q8WXI2
别名	CNKSR2;CNK2;KSR2;MAGUIN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn