

SNW1 Rabbit pAb

货号: B18892

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:200
理论分子量	61kDa
实测分子量	70kDa/75kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,293T,HeLa,HepG2,Mouse testis,Rat testis
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 300-536 of human SNW 1 (NP_036377.1).
序列	EAVEMRAQVERKMAQKEKEKHEEKLREMAQKARERRAGIKTHVEKEDGEARERDEIRHDRRKERQHDRNLSRAAPDKRS KLQRNENRDISIEVIALGVPNPRTSNEVQYDQRLFNQSKGMDSGFAGGEDEIYNVYDQAWRGGKDMAQSIYRPSKNLDK DMYGDDLEARIKTNRFPDKEFSGSDRRQRGREGPVQFEEDPFGLDKFLEEAKQHGGSKRPSDSSRPKEHEHEGKKRR KE

靶点信息

研究背景	This gene, a member of the SNW gene family, encodes a coactivator that enhances transcription from some Pol II promoters. This coactivator can bind to the ligand-binding domain of the vitamin D receptor and to retinoid receptors to enhance vitamin D-, retinoic acid-, estrogen-, and glucocorticoid-mediated gene expression. It can also function as a splicing factor by interacting with poly(A)-binding protein 2 to directly control the expression of muscle-specific genes at the transcriptional level. Finally, the protein may be involved in oncogenesis since it interacts with a region of SKI oncoproteins that is required for transforming activity. Alternative splicing results in multiple transcript variants.
基因ID	22938
基因名	SNW1
Swiss	Q13573
别名	SNW1;Bx42;NCOA-62;PRPF45;Prp45;SKIIP;SKIP;FUN20;SKIP1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn