

FAM13A Rabbit pAb

货号: B18897

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	77kDa/78kDa/80kDa/116kDa
实测分子量	116kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human FAM13 A (NP_055698.2).
序列	MGAGALAICQSKAAVRLKEDMKKIVAVPLNEQKDFTYQKLFGVSLQELERQGLTENGIPAVVWNIVEYLTQHGLTQEGLFRVNGNVKVVEQLRLKFESGVPVELGKGDVCSAASLLKFLRELPDSLITSALQPRFIQLFQDGRNDVQESSLRDLIKELPDTHYCLLKYL CQFLTKVAKHHVQNRMNVHNLATVFGPNC

靶点信息

研究背景	Predicted to enable GTPase activator activity. Predicted to be involved in regulation of small GTPase mediated signal transduction. Predicted to be located in cytosol. Implicated in chronic obstructive pulmonary disease.
基因ID	10144

基因名	FAM13A
Swiss	O94988
别名	FAM13A;ARHGAP48;FAM13A1

产品验证



Western blot analysis of FAM13A expressed in Jurkat,HeLa using FAM13A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn