

NGDN Rabbit pAb

货号: B18911

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	35kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,K-562,HepG2,Mouse testis,Mouse brain,Mouse spleen,Rat testis
细胞定位	Cell projection,Chromosome,Cytoplasm,Nucleus,axon,centromere,dendrite,filopodium,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 61-209 of human NGDN (NP_001036100.1).
序列	LMYLMDLTHLILDKASGGSLQGHDAVLRLVEIRTVLEKLRPLDQKLKYQIDKLIKTAVTGSLSENDPLRFKPHPSNMMSKLS SEDEEEDAEDDQSEASGKKS VKGVSKKYVPPRLVPVHYDETEAEREKKRLERAKRRALSSSVIREL

靶点信息

研究背景	
基因ID	25983
基因名	NGDN

Swiss	Q8NEJ9
别名	NGDN;C14orf120;CANu1;LCP5;NGD;lpd-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn