

# MYCT1 Rabbit pAb

货号: B18913

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 26kDa                                                                                                |
| 实测分子量 | 27kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Rat spleen,Rat lung                                                                                  |
| 细胞定位  | Nucleus                                                                                              |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                             |
|------|---------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 99-235 of human MYCT 1 (NP_079383.2).                         |
| 序列   | ASAPISQWSSSRRSRSSYTHGLNRTGFYRHSGCERRSNLSLASLTFQRQASLEQANSFPRKSSFRASTFHPFLQCPLPVE<br>TESQLVTLPSSNISPTISTSHSLSRPDYWSSNSLRVGLSTPPPPAYESIKAFPDS |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | May regulate certain MYC target genes, MYC seems to be a direct upstream transcriptional activator. Doe s not seem to significantly affect growth cell capacity. Overexpression seems to mediate many of the kno wn phenotypic features associated with MYC, including promotion of apoptosis, alteration of morphology, enhancement of anchorage-independent growth, tumorigenic conversion, promotion of genomic instabilit y, and inhibition of hematopoietic differentiation (By similarity. |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |            |
|-------|------------|
| 基因ID  | 80177      |
| 基因名   | MYCT1      |
| Swiss | Q8N699     |
| 别名    | MYCT1;MTLC |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)