

# SNRPD1 Rabbit pAb

货号: B18933

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IHC IF/ICC                                                                                        |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200                 |
| 理论分子量 | 13kDa                                                                                                |
| 实测分子量 | 13KDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | HeLa,293T,K-562,Mouse brain,Mouse spleen,Rat brain,Rat thymus                                        |
| 细胞定位  | Cytoplasm,Nucleus,cytosol                                                                            |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-96 of human SNRPD1 (NP_008869.1). |
| 序列   | MKLVRFLMKLSHETVTIELKNGTQVHGTITGVDVSMNTHLKAVKMTLKNREPQVQLETLSIRGNNIRYFILPDSLPLDTLLVDV<br>EPKVKSKKREAVA             |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                               |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a small nuclear ribonucleoprotein that belongs to the SNRNP core protein family. The protein may act as a charged protein scaffold to promote SNRNP assembly or strengthen SNRNP-SNRNP interactions through nonspecific electrostatic contacts with RNA. Two transcript variants encoding different isoforms have been found for this gene. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                 |
|-------|---------------------------------|
| 基因ID  | 6632                            |
| 基因名   | SNRPD1                          |
| Swiss | P62314                          |
| 别名    | SNRPD1;HsT2456;SMD1;SNRPD;Sm-D1 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)