

GOLPH4 Rabbit pAb

货号: B18943

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	81kDa
实测分子量	130kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,LO2
细胞定位	Endosome membrane,Golgi apparatus,Golgi stack membrane,Lipid-anchor,Membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 34-300 of human GOLPH4 (NP_055313.1).
序列	YELQTQLRKAEAVALKYQQHQESLSAQLQVVYEHRSRLEKSLQKERLEHKKAKEDFLVYKLEAQETLNKGRQDSNSRYSA LNVQHQM LKSQHEELKKQHSDLEEEHRKQGEDFSRTFNDHKQKYQLQQEKEQELSKLKETVYNLREENRQLRKAHQDI HTQLQDVKQQHKNNLLSEHEQLVVTLEDHKSALAAQAQTQVAEYKQLKDTLNRIPSLRKPDPAEQQNVTVQVAHSPQGYNTA REKPTREVQEVSRNNDVWQNHEAVPGRAE

靶点信息

研究背景	The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is a type II Golgi-resident protein. It may process proteins synthesized in the rough endoplasmic reticulum and assist in the transport of protein cargo through the Golgi apparatus.
------	--

基因ID	27333
基因名	GOLIM4
Swiss	O00461
别名	GOLIM4;GIMPC;GOLPH4;GPP130;P138

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn