

B3GNT5 Rabbit pAb

货号: B18957

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa
实测分子量	44kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,A-549
细胞定位	Golgi apparatus membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 36-378 of human B3GN T5 (NP_114436.1).
序列	DNHIVSHMKSYSRYRLINSYDFVNDTLSLKHTSAGPRYQYLINHKEKCQAQDVLLLLFVKTAPENYDRRSGIRRTWGNENY VRSQLNANIKTLFALGTPNPLEGEELQRKLAWEDQRYNDIIQQDFVDSFYNLTLKLLMQFSWANTYCPHAKFLMTADDDIF IHMPNLIEYLSLEQIGVQDFWIGRVHRGAPPIRDKSSKYVSYEMYQWPAYPDYTAGAAYVISGDVAAKVYEASQTLNSS LYIDDVFMGLCANKIGIVPQDHVFFSGEGKTPYHPCIYEKMMTSHGHLEDLQDLWKNATDPKVKTISKGFFGQIYCRLMKII LLCKISYVDTPCRAAFI

靶点信息

研究背景	This gene encodes a member of the beta-1,3-N-acetylglucosaminyltransferase family. This enzyme is a ty pe II membrane protein. It exhibits strong activity to transfer GlcNAc to glycolipid substrates and is identi fied as the most likely candidate for lactotriaosylceramide synthase. This enzyme is essential for the expr ession of Lewis X epitopes on glycolipids.
------	---

基因ID	84002
基因名	B3GNT5
Swiss	Q9BYG0
别名	B3GNT5;B3GN-T5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn