

TSPAN15 Rabbit pAb

货号: B18970

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa
实测分子量	33kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse spleen
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-220 of human TSPA N15 (NP_036471.1).
序列	GVVALTFRNQITIDFLNDNIRRGIIENYYDDLDFKNIMDFVQKKFKCCGGEDYRDWSKNQYHDCSAPGPLACGVPTCCIRN TTEVVNTMCGYKTIDKERFSVQDVIYVRGCT

靶点信息

研究背景	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the te traspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regul ation of cell development, activation, growth and motility. The use of alternate polyadenylation sites has been found for this gene.
基因ID	23555

基因名	TSPAN15
Swiss	O95858
别名	TSPAN15;2700063A19Rik;NET-7;NET7;TM4SF15

产品验证



Western blot analysis of TSPAN15 expressed in Mouse spleen using TSPAN15 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn