

ZNF703 Rabbit pAb

货号: B18989

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	58kDa
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 160-290 of human ZNF703 (NP_079345.1).
序列	GSGGGDSRKDSGSSSVSSTSSSSSSSPGDKAGFRVPSAACPPFPHPGAPVSASSSSSSPGGSRGGSPHHSDCKNGGGVGGGELDKKDQEPKPSPEPAAVSRGGGGEPGAHGGAESGASGRKSEPPSALVGA

靶点信息

研究背景	Transcriptional corepressor which does not bind directly to DNA and may regulate transcription through recruitment of histone deacetylases to gene promoters. Regulates cell adhesion, migration and proliferation. May be required for segmental gene expression during hindbrain development.
基因ID	80139

基因名	ZNF703
Swiss	Q9H7S9
别名	ZNF703;NLZ1;ZEPP01;ZNF503L;ZPO1

产品验证



Western blot analysis of ZNF703 expressed in MCF7 using ZNF703 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn