

ILKAP Rabbit pAb

货号: B19004

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	42kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,Mouse spleen,Mouse testis
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human ILKAP (NP_110395.1).
序列	MDLFGDLPEPERSPRPAAGKEAQKGPLLFDLPPASSTDGSGGGPLLFDLPPASSGDSGSLATSISQMVKTEGKGAKRK TSEEEKNGSEELVEKKVCKASSVIFGLKGYVAERKGEREEMQDAHVILNDITECRPPSSLITRVSYFAVFDGHGGIRASKFA AQN LHQNLIRKFPKGDVISVEKTVKRCLLDTFKHTDEEFLKQASSQKPAWKDGSTATCVLAVDNILYIANLGDSRAILCRYN EESQKHAALSLSKEHNPTQYEERMRIQAGGNVRDGRVLGVLEVSRSIGDGQYKRCGVTSVPDIRRCQLTPNDRFILLAC DGLFKVFTPEEAVNFILSCLEDEKIQTREGKSAADARYEACNRLANKAVQRGSADNVTVMVVRIGH

靶点信息

研究背景	The protein encoded by this gene is a protein serine/threonine phosphatase of the PP2C family. This prote in can interact with integrin-linked kinase (ILK/ILK1), a regulator of integrin mediated signaling, and regul ate the kinase activity of ILK. Through the interaction with ILK, this protein may selectively affect the sign aling process of ILK-mediated glycogen synthase kinase 3 beta (GSK3beta), and thus participate in Wnt si gnaling pathway.
------	--

基因ID	80895
基因名	ILKAP
Swiss	Q9H0C8
别名	ILKAP;ILKAP2;ILKAP3;PP2C-DELTA;PPM1O

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn