

IGSF8 Rabbit pAb

货号: B19007

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	40kDa/55kDa/65kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji
细胞定位	Cell membrane,Single-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 28-350 of human IGSF8 (NP_443100.1).
序列	REVLVPEGPLYRVAGTAVSISCNVTGYEGPAQQNFEWFLYRPEAPDTALGIVSTKDTQFSYAVFKSRVVAGEVQVQRLQG DAVVLKIARLQAQDAGIYECHTPSTDTRYLGSYSGKVELRVLDPVLQVSAAPPGPRGRQAPTSPPRMTVHEGQELALGCLA RTSTQKHTHLAVSFGRSVPEAPVGRSTLQEYVVGIRSDLAVEAGAPYAERLAAGELRLGKEGTDYRMYVVGGAQAGDAGT YHCTAAEWIQDPDGSWAQIAEKRAVLAHVDVQTLSSQLAVTVGPGERRIGPGEPLLELLCNVSGALPPAGRHAAYSVGWE MAPA

靶点信息

研究背景	This gene encodes a member the EWI subfamily of the immunoglobulin protein superfamily. Members of this family contain a single transmembrane domain, an EWI (Glu-Trp-Ile)-motif and a variable number of immunoglobulin domains. This protein interacts with the tetraspanins CD81 and CD9 and may regulate their role in certain cellular functions including cell migration and viral infection. The encoded protein may also function as a tumor suppressor by inhibiting the proliferation of certain cancers. Alternate splicing results in multiple transcript variants that encode the same protein.
基因ID	93185
基因名	IGSF8
Swiss	Q969P0
别名	IGSF8;CD316;CD81P3;EWI-2;EWI2;KCT-4;LIR-D1;PGRL

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn