

IQGAP3 Rabbit pAb

货号: B19050

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	184kDa
实测分子量	170kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	cell-cell junction,cytosol,lateral plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1437-1631 of human IQ GAP3 (NP_839943.2).
序列	QRRVLRNLRRLREALGLVSARNGYQGLVDELAKDIRNQHRHRHRRKAELVKLQATLQGLSTKTTFYEEQGDYYSQYIRACLD HLAPDSKSSGKGGKKQPSLHYTAAQLLEKGVLEIEDLPASHFRNVIFDITPGDEAGKFEVNAKFLGVDMERFQLHYQDLLQ LQYEGVAVMKLFNKAKVNVNLLIFLLNKKFLRK

靶点信息

研究背景	Enables calmodulin binding activity and myosin VI light chain binding activity. Predicted to be involved in regulation of actin cytoskeleton organization. Predicted to act upstream of or within several processes, including intracellular signal transduction; positive regulation of macromolecule metabolic process; and positive regulation of mammary gland epithelial cell proliferation. Predicted to be located in cytosol. Predicted to be active in cytoplasm.
------	--

基因ID	128239
基因名	IQGAP3
Swiss	Q86VI3
别名	IQGAP3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn