

GAS7 Rabbit pAb

货号: B19054

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	39kDa/47kDa/54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A375,Mouse brain,Mouse testis,Rat brain
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-412 of human GAS7 (NP_001124303.1).
序列	MVPPPPGEESQTVILPPGWQSYLSPQGRRYYVNTTTNETTWERPSSSPGIPASPGSHRSSLPPTVNGYHASGTPAHPPETA HMSVRKSTGDSQNLGSSSPSKKQSKENTITINCVTFFHPDTMPEQQLLKPTESYCDYFWADKKDPQGNGTVAGFELL QKQLKGQMCKEMSEFIRERIKIEEDYAKNLAKLSQNSLASQEEGSLGEAWAQVKKSLADEAEVHLKFSAKLHSEVEKPL MNFRENFKKDMKKCDHHIADLRKQLASRYASVEKARKALTERQRDLEMKTQQLEIKLSNKTEEDIKKARRKSTQAGDDL MRCVDLYNQAQSKWFEEMVTTTLELERLEVERVEMIRQHLCQYTLRHETDMFNQSTVEPVDQLLRKVDPKADRELWVRE HKTGNIRPVDMEI

靶点信息

研究背景	Growth arrest-specific 7 is expressed primarily in terminally differentiated brain cells and predominantly in mature cerebellar Purkinje neurons. GAS7 plays a putative role in neuronal development. Several transcript variants encoding proteins which vary in the N-terminus have been described.
------	---

基因ID	8522
基因名	GAS7
Swiss	O60861
别名	GAS7;MLL/GAS7

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn