

# RNF31 Rabbit pAb

货号: B19064

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                            |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB IF/ICC                                                                                              |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:100                                               |
| 理论分子量 | 57kDa/102kDa/119kDa                                                                                    |
| 实测分子量 | 120kDa                                                                                                 |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | BT-474,PC-3,Mouse brain                                                                                |
| 细胞定位  | Cytoplasm                                                                                              |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                                                                                                   |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human RNF31 (NP_060469.4).                                                                                               |
| 序列   | AAGACPEEIFSALQYSGTEVPLQWLRSELPYVLEMVAELAGQQDPGLGAFSCQEARRAWLDRHGNNLDEAVEECVTRRRK<br>VQELQSLGFGPEEGSLQALFQHGGDVSRALTELQRQRLEPFRQRLWDSGPEPTPSWDGPDKQSLVRRLLAVYALPSW<br>GRAELALSLLQETPRNYELGDVVEAVRHSQDRAFLRRLLAQECA |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                       |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. The encoded protein is the E3 ubiquitin-protein ligase component of the linear ubiquitin chain assembly complex. Two transcript variants encoding different isoforms have been found for this gene. |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                       |
|-------|-----------------------|
| 基因ID  | 55072                 |
| 基因名   | RNF31                 |
| Swiss | Q96EP0                |
| 别名    | RNF31;HOIP;Paul;ZIBRA |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)