

NFIL3 Rabbit pAb

货号: B19096

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:100
理论分子量	51kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 133-462 of human NFIL 3 (NP_005375.2).
序列	AYAQEIQKLSNSTAVYFQDYQTSKSNVSSFVDEHEPSMVSSSCISVIKHSPQSSLSDVSEVSSVEHTQESSVQGSCRSPEN KFQIIKQEPMELESYTREPRDDRGSYTASIYQNYMGNSFSGYSHSPPLLQVNRSSSNSPRTSETDDGVVGKSSDGEDEQQ VPKGPISPVELKHVHATVVKVPEVNSSALPHKLRIKAKAMQIKVEAFDNEFEATQKLSSPIDMTSKRHFLEKHSAPSMVH SSLTPFSVQVTNIQDWSLKSEHWHQKELSGKTQNSFKTGVVEMKDSGYKVSDPENLYLKQGIANLSAEVVSLKRLIATQPI SASDSG

靶点信息

研究背景	The protein encoded by this gene is a transcriptional regulator that binds as a homodimer to activating tr anscription factor (ATF) sites in many cellular and viral promoters. The encoded protein represses PER1 a nd PER2 expression and therefore plays a role in the regulation of circadian rhythm. Three transcript varia nts encoding the same protein have been found for this gene.
------	---

基因ID	4783
基因名	NFIL3
Swiss	Q16649
别名	NFIL3;E4BP4;IL3BP1;NF-IL3A;NFIL3A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn