

MESDC2 Rabbit pAb

货号: B19127

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HepG2,SGC-7901,Mouse testis
细胞定位	Endoplasmic reticulum
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human MESDC 2 (NP_055969.1).
序列	MAASRWARKAVVLLCASDLLLLLLLLPPPGSCAAEGSPGTPDESTPPPRKKKKDIRDYNDADMARLLEQWEKDDDDIEEGD LPEHKRPSAPVDFSKIDPSKPESILKMTKKGKTLMMFVTVSGSPTEKETEEITSLWQGSLFNANYDVQRFIVGSDRAIFMLR DGSYaweikDFLVGQDRCADVTLEGQVYPGKGGGSKEKNKTKQDKGKKKKEGDLKSRSSKEENRAGNKREDL

靶点信息

研究背景	Chaperone specifically assisting the folding of beta-propeller/EGF modules within the family of low-density lipoprotein receptors (LDLRs). Acts as a modulator of the Wnt pathway through chaperoning the coreceptors of the canonical Wnt pathway, LRP5 and LRP6, to the plasma membrane. Essential for specification of embryonic polarity and mesoderm induction. Plays an essential role in neuromuscular junction (NMJ) formation by promoting cell-surface expression of LRP4 (By similarity). May regulate phagocytosis of apoptotic retinal pigment epithelium (RPE) cells (By similarity).
------	---

基因ID	23184
基因名	MESDC2
Swiss	Q14696
别名	MESDC2;BOCA;MESD

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn