

CPT1C Rabbit pAb

货号: B19163

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	81kDa/89kDa/90kDa
实测分子量	91kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse testis,Mouse kidney,Rat brain,Rat testis,Rat kidney,L929
细胞定位	Cell junction,Cell projection,Endoplasmic reticulum,Mitochondrion outer membrane,Multi-pass membrane protein,axon,dendrite,synapse
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 270-400 of human CPT1C (NP_001129524.1).
序列	LYRHLNRQEIPPTLLMGMRPLCSAQYEKIFNTTRIPGVQKDYIRHLHDSQHVAVFHRGRFFRMGTHSRNSLLSPRALEQQFQRILDDPSPACPHEEHLAALTAAPRGTTWAQVRTSLKTQAAEALEAVEGA

靶点信息

研究背景	This gene encodes a member of the carnitine/choline acetyltransferase family. The encoded protein regulates the beta-oxidation and transport of long-chain fatty acids into mitochondria, and may play a role in the regulation of feeding behavior and whole-body energy homeostasis. Alternatively spliced transcript variants encoding multiple protein isoforms have been observed for this gene.
------	---

基因ID	126129
基因名	CPT1C
Swiss	Q8TCG5
别名	CPT1C;CATL1;CPT1-B;CPT1P;CPTI-B;CPTIC;SPG73

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn