

# DLGAP2 Rabbit pAb

货号: B19181

产品信息

|       |                                                                                                                 |
|-------|-----------------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                           |
| 宿主    | Rabbit                                                                                                          |
| 克隆性   | Polyclonal                                                                                                      |
| 预测反应  |                                                                                                                 |
| 应用    | WB                                                                                                              |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                                       |
| 理论分子量 | 108kDa/115kDa/117kDa                                                                                            |
| 实测分子量 | 117kDa                                                                                                          |
| 形式    | Liquid                                                                                                          |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.            |
| 偶联物   | Unconjugated                                                                                                    |
| 阳性对照  | U-87MG,293T                                                                                                     |
| 细胞定位  | Cell junction,Cell membrane,Peripheral membrane protein,postsynaptic cell membrane,postsynaptic density,synapse |
| 纯化    | Affinity purification                                                                                           |

抗原信息


|      |                                                                                                                                                                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 550-800 of human DLGAP2 (NP_004736.2).                                                                                                                                                 |
| 序列   | TAQSSTESTQDAYQDSRAQRMSPWPQDSRGLYNSTDSLDSNKAMNLALETAAQRLPESQSSSVRTSDKAILVSKAE<br>ELLKSRCSSIGIQDSEFPEHQPYPRSDVETATDSDTESRGLREYHSGVGQVEDEKRHGRFKRSNSVTAAVQADLELEGFP<br>GHITTEDKGLQFGSSFQRHSEPSTPTQYSAVRTVRTQGLFSYREDYRTQVDTSTLPPDPWLEPAIDTVETGRMSPCRRD<br>GSWFLKLLHAETK |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                       |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The product of this gene is a membrane-associated protein that may play a role in synapse organization and signalling in neuronal cells. This gene is biallelically expressed in the brain, however, only the paternal allele is expressed in the testis (PMID:18055845). Alternatively spliced transcript variants encoding different isoforms have been identified. |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                       |
|-------|---------------------------------------|
| 基因ID  | 9228                                  |
| 基因名   |                                       |
| Swiss | Q9P1A6                                |
| 别名    | DLGAP2;C8orf68;DAP2;ERICH1-AS1;SAPAP2 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)