

# FCRL5 Rabbit pAb

货号: B19185

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 14kDa/17kDa/65kDa/83kDa/106kDa                                                                       |
| 实测分子量 | 106kDa                                                                                               |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Jurkat,PC-12,SP2/0                                                                                   |
| 细胞定位  | Cell membrane,Single-pass type I membrane protein                                                    |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                          |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 16-190 of human FCRL5 (NP_001182317.1).                                                                    |
| 序列   | QFARTPRPIIFLQPPWTTVFQGERVTLTCKGFRFYSPQKTKWYHRYLGKEILRETPDNILEVQESGEYRCQAQGSPLSSPVH<br>LDFSSASLILQAPLSVFEGDSVVLRCRAKAEVTLNNTIYKNDNVLAFLNKRRTDFHIPHACLKDNGAYRCTGYKESCCPVSSN<br>TVKIQVQEPFT |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a member of the immunoglobulin receptor superfamily and the Fc-receptor like family. This gene and several other Fc receptor-like gene members are clustered on the long arm of chromosome 1. The encoded protein is a single-pass type I membrane protein and contains 8 immunoglobulin-like C2-t ype domains. This gene is implicated in B cell development and lymphomagenesis. Alternatively spliced t ranscript variants encoding different isoforms have been identified. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                                              |
|-------|----------------------------------------------|
| 基因ID  | 83416                                        |
| 基因名   | FCRL5                                        |
| Swiss | Q96RD9                                       |
| 别名    | FCRL5;BXMAS1;CD307;CD307e;FCRH5;IRTA2;PRO820 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)