

# KLF17 Rabbit pAb

货号: B19204

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 | 42kDa                                                                                                |
| 实测分子量 | 43kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | U-87MG,HepG2,MCF7,Mouse brain                                                                        |
| 细胞定位  | Nucleus                                                                                              |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                             |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 210-389 of human KLF17 (NP_775755.3).                                                                         |
| 序列   | QMLPPQDAHDLGMPPAESQSLVLVLSQDSLVSQPDSQEGPFLPEQPGPAPQTVEKNSRPQEGTGRRGSSEARPYCCNY<br>ENCGKAYTKRSHLVSHQRKHTGERPYSCNWESCSWSFFRSDELRRHMRVHTRYRPYKCDQCSREFMRSDDLKQHQKT<br>HRPGPSDPQANNNGEQDSPPAAGP |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Transcription repressor that binds to the promoter of target genes and prevents their expression. Acts as a negative regulator of epithelial-mesenchymal transition and metastasis in breast cancer. Specifically binds the 5'-CACCC-3' sequence in the promoter of ID1, a key metastasis regulator in breast cancer, and repress its expression. May be a germ cell-specific transcription factor that plays important roles in spermatid differentiation and oocyte development (By similarity). |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                     |
|-------|---------------------|
| 基因ID  | 128209              |
| 基因名   | KLF17               |
| Swiss | Q5JT82              |
| 别名    | KLF17;ZNF393;Zfp393 |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)