

CABP7 Rabbit pAb

货号: B19223

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	24kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse lung,Mouse kidney,Rat lung,Rat kidney
细胞定位	Cell membrane,Cytoplasm,Golgi apparatus,Single-pass type IV membrane protein,perinuclear region,trans-Golgi network membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human CABP 7 (NP_872333.1).
序列	PKLSTSGIPEKFHGTDFDTVFWKCDMQKLTVDELKRLLYDTFCEHLSMKDIENIIMTEEEESHLGTAEECPVDVETCSNQQI

靶点信息

研究背景	Negatively regulates Golgi-to-plasma membrane trafficking by interacting with PI4KB and inhibiting its activity.
基因ID	164633

基因名	CABP7
Swiss	Q86V35
别名	CABP7;CALN2

产品验证



Western blot analysis of CABP7 expressed in Mouse brain,Rat lung,Mouse lung using CABP7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn