

RND1 Rabbit pAb

货号: B19224

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Rat brain
细胞定位	Cell membrane,Cytoplasm,Cytoplasmic side,Lipid-anchor,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-232 of human RND1 (NP_055285.1).
序列	EQGCAIAKQLGAEIYLEGSAFTSEKSIHSIFRTASMLCLNKPSLPQKSPVRSLSKRLLHLP SRSELISSTFKKEKAKSCSIM

靶点信息

研究背景	This gene encodes a protein that belongs to the Rho GTPase family. Members of this family regulate the organization of the actin cytoskeleton in response to extracellular growth factors. A similar protein in rat interacts with a microtubule regulator to control axon extension.
基因ID	27289

基因名	RND1
Swiss	Q92730
别名	RND1;ARHS;RHO6;RHOS

产品验证



Western blot analysis of RND1 expressed in Mouse brain,Mouse liver,Rat brain using RND1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn