

PFN4 Rabbit pAb

货号: B19229

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	14kDa
实测分子量	14kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-129 of human PFN4 (NP_955378.1).
序列	MSHLQSLLLDTLLGTHVDSAALIKIQERSLCVASPGFNVTPSDVRTLVNGFAKNPLQARREGLYFKGKDYRCVRADEYSLYAKNENTGVVVVKTHLYLLVATYTEGMYPsicVEATESLGDYLRKKGS

靶点信息

研究背景	Binds to phosphatidylinositol 3-phosphate (PtdIns(3P, phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5P2 , phosphatidylinositol 4-phosphate (PtdIns(4P and phosphatidic acid (PA. Does not bind to actin, contrary to other family members.
基因ID	375189

基因名	PFN4
Swiss	Q8NHR9
别名	PFN4

产品验证



Western blot analysis of PFN4 expressed in Mouse heart,Rat heart using PFN4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn